

Myelo-Nate[®]

DEHP
Free

Neonatal/Pediatric CSF Sampling Kit

...minimizes trauma from cerebrospinal space access during neonatal/pediatric lumbar punctures.



The needle's crystal clear hub provides quick visual detection of fluid flow.

The short, bevel design with Quincke-type point spreads and displaces the dura mater.

This needle reduces overall penetration forces, yet still imparts a positive feel when penetrating the dura.

UTAH MEDICAL
PRODUCTS INC.



- ◆ Choice of 22 ga. or 25 ga. lumbar puncture needle with stylet
- ◆ Thin wall cannula for better flow
- ◆ Short bevel point geometry
- ◆ Integrally ground stylet
- ◆ Contoured hubs for better tactile sensitivity
- ◆ Crystal clear view hub
- ◆ Color coded stylet cap
- ◆ Marker for bevel-up orientation

Product Description

The Myelo-Nate Kit is designed to provide all the essentials necessary for performing a lumbar puncture procedure for diagnostic testing of spinal fluid on neonatal and pediatric patients.

This procedure removes cerebral spinal fluid (CSF), which is analyzed for the presence of bacteria, glucose, protein, and blood cells.

Complete Kit Contains:

- One Myelo-Nate LP needle with stylet
- One Point Lok needle retention device
- Four graduated specimen vials
- Povidone-Iodine swabsticks
- Gauze pads
- Towel
- Fenestrated drape
- Spot bandage

Myelo-Nate Needle

<u>Catalog Number</u>	<u>Myelo-Nate Description</u>	<u>Qty per Box</u>
4011019	22 ga x 1" Needle	30 per box
4011529	22 ga x 1.5" Needle	30 per box
4011519	25 ga x 1" Needle	30 per box

Myelo-Nate Complete Kit

<u>Catalog Number</u>	<u>Myelo-Nate Description</u>	<u>Qty per Box</u>
4011015	22 ga x 1" Kit	10 per box
4011525	22 ga x 1.5" Kit	10 per box
4011515	25 ga x 1" Kit	10 per box

Additional Advantages:

- The complete tray provides all the necessary items for the sterile collection of CSF.
- The LP Needle is packaged separately in 1 inch or 1.5 inch lengths.
- The thin wall cannula design enhances flow of CSF.
- The crystal clear hub allows early identification of fluid flow.
- The needle has an integrally ground stylet and a marker for bevel up orientation.

UTAH MEDICAL
PRODUCTS INC.



Gesco™

UNITED STATES
7043 South 300 West
Midvale, Utah 84047
(801) 566-1200
Toll Free 1-800-533-4984
Fax (801) 566-2062
www.utahmed.com

 EUROPE
Utah Medical Products, Ltd.
Athlone Business & Technology Park
Athlone, County Westmeath, Ireland
Tel: +353-90-647-3932

P/N 58283 rev 063008