

PRESORT STANDARD
U. S. POSTAGE
PAID
PERMIT NO. 3280
SALT LAKE CITY, UT

**UTAH MEDICAL
PRODUCTS INC.**
7043 South 300 West ♦ Midvale, Utah 84047
800-533-4984 ♦ 801-566-1200
www.utahmed.com


An infection is the most unwanted...



...when he is already fighting for his life.

Deltran[®] reduces neonatal sepsis.

Nosocomial infections can be fatal to your tiniest patients. That is why Utah Medical offers a solution to help prevent pathogen entry to the bloodstream: The Deltran[®] closed system.



- Closed system reduces nosocomial infections and line sepsis by preventing pathogen entry to the bloodstream.
- Blunt sampling cannula eliminates needle-stick exposure.
- Clearing volume tubing provides a saline-free sample for lab analysis.
- Reservoir syringe preserves precious blood volume for returning to the neonate.
- Integral dual function Snap-Tab allows dynamic response testing.

For more information or to evaluate Deltran[®] in your facility call
1-800-533-4984 or visit
www.utahmed.com/NICU.htm

